



Excellence in Ultrasonic Cleaning

A Message from Management

2005 at Soniclean is shaping up to be a very eventful year. We will be attending many exhibitions nationally and internationally.

A major milestone for our team at Soniclean is the launch of our Irrigated Ultrasonic Cleaner with Machine Reporting System (MRS). This technology enables the hospitals to record vital information of the cleaning process for instruments used on patients, back to the instrument cleaning and sterilization process. (See Feature Story below)

Our Sales and Marketing team will be working hard to gain broader representation for Soniclean's range of Ultrasonic Cleaners in the international market. Some of our target

markets are Asia, South Africa and with our CE mark compliance and the appointment of Rod Ruston of Priory Analysts as our Authorised Representative we are well on our way to Europe & UK!

We listen to our customers and value all forms of feedback. As a direct result of our customer's requirements and feedback we are continually looking at new ideas and have a number of new products and improvements coming through our R&D pipeline.

For detailed information on our Soniclean product range or for further information, please visit our web site www.soniclean.com.au

To request information, obtain quotes or provide feedback, email us at sales@soniclean.com.au

Upcoming Exhibitions:

- Decontamination Sciences Congress 2005 London. 30 March – 1 April
- South Australian SRACA Conference & Trade Show, Adelaide 29-30 April
- OTA Annual Conference & ACORN AGM, Canberra 14-16 April

Feature story:

- Irrigated Ultrasonic Cleaner with Machine Reporting System (MRS)

Inside this issue:

- Staff Profile – Sales & Export
Aaron Ng & Jodie Sabato
- Soniclean Delivers first Pre-Sterilizer with Machine Reporting System
- 31st Australian Dental Congress & Exhibition
- Soniclean Provides Contract Cleaning Service

Soniclean New Product Launch - Irrigated Ultrasonic Cleaner with Machine Reporting System (MRS)

Soniclean's Innovative Irrigated Ultrasonic Cleaner with Machine Reporting System (MRS) enables the hospitals to record vital information of the cleaning process for instruments used on a patient, back to the instrument cleaning and sterilization process.

Cannulated instruments are difficult and time consuming to clean. Using a Soniclean Irrigator is the ideal way to reduce the risk of infection and the time your staff takes to manually clean. The irrigator cleans and decontaminates all types of hollow instruments prior to sterilization. Ultrasonics have been shown to be much more effective than manual cleaning. It is common knowledge that surgical instruments must be clean of debris before sterilization can be effective.

Soniclean's Ultrasonic Irrigators use Soniclean's unique patented *Pulse Swept Power* technology to critically clean with superior penetration,



precision and reliability. This technology uses pulsed ultrasonics to ensure powerful, reliable, uniform cleaning of internal lumens and fine parts without the risk of damage to delicate instruments. The pulsed irrigation of clean fresh ultrasonically charged fluid produces maximum turbulence for cleaning the inner lumen. This combination of pulsed ultrasonic cleaning and pulsed flushing of the lumened instrument with clean fresh fluid means your instruments are cleaned faster and more thoroughly.

An instrument cannot be completely clean if it is flushed with recycled fluid as is the case with some other products on the market.

Soniclean's MRS logs each clean process internally; data is gathered and processed electronically during the clean process. Here are just a few of the many excellent advantages of Soniclean's MRS: Data is easy to file away as it is stored in computer file format; Data can be optionally printed out, the MRS does not insist on printed reports so it is also environmentally friendly on resources; The MRS can be accessed by either a PC connected directly to the cleaner or the cleaner can be permanently connected to a Local Area Network (LAN) and accessed by any networked client PC's; Provides a method of validating the machine against the original Soniclean validation Certificate carried out at the factory.

Please contact Soniclean for more information on MRS

31st Australian Dental Exhibition

The 31st Australian Dental Congress & Exhibition was held at the Adelaide Convention Centre, Adelaide, Australia on the 3rd - 7th of March 2005.


Soniclean was very well represented at the Exhibition and would like to thank our distributors Gunz Dental, Independent Dental Supplies, Davis Dental Depot and Siltex Australia who displayed our Soniclean Benchtop range of Ultrasonic Cleaners.



Aaron & Jodie with Christine from Gunz Dental



Aaron & Jodie with the team from Independent Dental Supplies



Julian, Jodie & Aaron with Jack Davis from Davis Dental Depot



Aaron & Jodie with the team from Siltex Australia

Soniclean Provides Contract Cleaning Service

Soniclean provides a contract ultrasonic cleaning service to our customers. We can help you by:

- ✚ Cleaning your articles for you,
- ✚ Providing advice to help you solve your cleaning problems,
- ✚ Renting immersible ultrasonic cleaning machines to you for your cleaning applications.

Please contact us for information on this excellent service

Staff Profile – Sales & Export Aaron Ng & Jodie Sabato

Aaron Ng has been Soniclean's Sales & Export Manager for the past 15 years. In this time he has built strong distributor networks in Australia and in over 60 countries worldwide such as Asia, Middle East, Far East, South Africa, Europe and North America.

Aaron's knowledge and expertise has enabled him to develop a very successful web based sales and marketing strategy, which has allowed Soniclean to grow at a rapid rate and to become very competitive in the international world of Ultrasonic Cleaners.



Jodie Sabato, who joined Soniclean in July 2004 as Sales & Export Assistant, supports Aaron. Jodie has had 5 years experience in export and sales. In 2003 she successfully completed the Australian Institute of Export's 'Certificate in International Trade Operations' course.

If you are interested in distributing our products please do not hesitate to contact Aaron by email or telephone aaron@soniclean.com.au, mobile 0403067181 or Jodie at jodie@soniclean.com.au, mobile 0433119308 or alternatively on the Soniclean office telephone +61 8 8234 8398.



Regions Soniclean currently has distributors in

Soniclean Delivers 1st Pre-Sterilizer with Machine Reporting System

We have delivered our first Soniclean Irrigated Pre-Sterilizer Machine with Lifter unit using the Soniclean Machine Reporting System (MRS). It has been installed at the Repatriation General Hospital in South Australia.

The Hospital's Central Sterile Supply Department (CSSD) will now be able to record and download valuable information for each cleaning cycle of the machine. The MRS will report status of recommended and pre tuned performance parameters and any irregularities with its unique clean number with time and date and much more information.

Soniclean's MRS is available as an optional feature for the Soniclean range of Ultrasonic Irrigators and Pre Sterilizers.

The Women's and Children's Hospital in Adelaide have ordered an Ultrasonic Irrigator with the MRS for precise critical cleaning of their cannulated instruments.



TGA Registration
Class 1 Certification
ARTG # 104498



Head Office


Soniclean Pty Ltd
38 Anderson Street
Thebarton, SA, 5031
Australia
Ph: + 61 8 8234 8398
Fax: + 61 8 8234 8391
ABN: 89 057 730 917

Malaysia Office

Sonic Solutions Sdn Bhd
13 Lorong 1/1
Taman Bukit Wuri Indah
70300 Seremban
N.S. Malaysia
Ph: + 60 6 631 6133
Fax: + 60 6 632 7507

Further Information

Email:
sales@soniclean.com.au
Web:
www.soniclean.com.au



AS/NZS ISO 9001:2000
Approved by BVQI