

Soniclean Ultrasonic Irrigators and Washers

Preventative Maintenance Contract Inclusions

1. All items listed under “General Service Inclusions”
2. Validation of ultrasonic performance
3. Soniclean Care – a discount on the monthly supply of detergent
4. Option to select number of visits per year
5. Onsite performance qualification (PQ)
6. Detailed PQ report including cleaning efficacy with certificate
7. Immediate access to the latest software and hardware upgrades
8. Free loan unit (freight charges apply)

General Service Inclusions

1. Unit inspection including cabinet, lid, tank, and fascia
2. Inspection of the display terminal, plumbing system, control module, water level sensor, venturi soap injector, lifter mechanism, and external connections
3. Inspection of printer and LAN connections
4. Check built-in diagnostics reports
5. Annual electrical safety check
6. Annual re-calibration

